DRI OCT Série Triton

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- “ Ikuno Y, Kawaguchi K, Nouchi T, Yasuno Y., “Choroidal thickness in healthy Japanese subjects.”, *Invest Ophthalmol Vis Sci*. 2010 Apr;51(4):2173-6.
- “ Hirata M, Tsujikawa A, Matsumoto A, Hangai M, Ooto S, Yamashiro K, Akiba M, Yoshimura N., “Macular choroidal thickness and volume in normal subjects measured by swept-source optical coherence tomography.”, *Invest Ophthalmol Vis Sci*. 2011 Jul 1;52(8):4971-8.
- “ Ikuno Y, Maruko I, Yasuno Y, Miura M, Sekiryu T, Nishida K, Iida T., “ Reproducibility of retinal and choroidal thickness measurements in enhanced depth imaging and high-penetration optical coherence tomography.”, *Invest Ophthalmol Vis Sci*. 2011 Jul 25;52(8):5536-40.
- “ Maruko I, Iida T, Sugano Y, Oyamada H, Sekiryu T., “Morphologic choroidal and scleral changes at the macula in tilted disc syndrome with staphyloma using optical coherence tomography.”, *Invest Ophthalmol Vis Sci*. 2011 Nov 11;52(12):8763-8.
- “ Ohno-Matsui K, Akiba M, Moriyama M, Ishibashi T, Tokoro T, Spaide RF., “Imaging retrobulbar subarachnoid space around optic nerve by swept-source optical coherence tomography in eyes with pathologic myopia.”, *Invest Ophthalmol Vis Sci*. 2011 Dec 28;52(13):9644-50.
- “ Usui S, Ikuno Y, Miki A, Matsushita K, Yasuno Y, Nishida K., “Evaluation of the choroidal thickness using high-penetration optical coherence tomography with long wavelength in highly myopic normal-tension glaucoma.”, *Am J Ophthalmol*. 2012 Jan;153(1):10-6.e1.
- “ Tsuchiya K, Moriyama M, Akiba M, Tamura Y, Ohno-Matsui K., “Development of peripapillary venous loop in an eye with a small optic disc.”, *Int Ophthalmol*. 2012 Apr;32(2):171-5. Epub 2012 Feb 10.
- “ Spaide RF, Akiba M, Ohno-Matsui K., “Evaluation of peripapillary intrachoroidal cavitation with swept source and enhanced depth imaging optical coherence tomography.”, *Retina*. 2012 Jun;32(6):1037-44.
- “ Ohno-Matsui K, Akiba M, Moriyama M, Shimada N, Ishibashi T, Tokoro T, Spaide RF., “Acquired Optic Nerve and Peripapillary Pits in Pathologic Myopia.”, *Ophthalmology* 2012; 119: 1685-1692
- “ Jirattanasopa P, Ooto S, Tsujikawa A, Yamashiro K, Hangai M, Hirata M, Matsumoto A, Yoshimura N., “Assessment of Macular Choroidal Thickness by Optical Coherence Tomography and Angiographic Changes in Central Serous Chorioretinopathy.”, *Ophthalmology* 2012; 119: 1685-1692
- “ Usui S, Ikuno Y, Akiba M, Maruko I, Sekiryu T, Nishida K, Iida T., “Circadian changes in subfoveal choroidal thickness and the relationship with circulatory factors in healthy subjects.”, *Invest Ophthalmol Vis Sci*. 2012 Apr 24;53(4):2300-7.
- “ Ohno-Matsui K, Akiba M, Moriyama M, Ishibashi T, Hirakata A, Tokoro T., “Intrachoroidal Cavitation in Macular Area of Eyes With Pathologic Myopia.”, *Am J Ophthalmology* 2012; 154:382-393
- “ Ellabban AA, Tsujikawa A, Matsumoto A, Ogino K, Hangai M, Ooto S, Yamashiro K, Akiba M, Yoshimura N., “Macular choroidal thickness and volume in eyes with angioid streaks measured by swept source optical coherence tomography.”, *Am J Ophthalmol*. 2012 Jun;153(6):1133-1143.e1.
- “ Maruko I, Iida T, Sugano Y, Oyamada H, Akiba M, Sekiryu T., “Morphologic analysis in pathologic myopia using high-penetration optical coherence tomography.”, *Invest Ophthalmol Vis Sci*. 2012 Jun 20;53(7):3834-8
- “ Ohno-Matsui K, Akiba M, Modegi T, Tomita M, Ishibashi T, Tokoro T, Moriyama M., “Association between Shape of Sclera and Myopic Retinochoroidal Lesions in Patients with Pathologic Myopia.”, *Invest Ophthalmol Vis Sci*. 2012 Sep 7;53(10):6046-61
- “ Ohno-Matsui K, Akiba M, Ishibashi T, Moriyama M., “Observations of Vascular Structures within and Posterior to Sclera in Eyes with Pathologic Myopia by Swept-Source Optical Coherence Tomography.”, *Invest Ophthalmol Vis Sci*. 2012 Oct 19;53(11):7290-8
- “ Ellabban AA, Tsujikawa A, Matsumoto A, Yamashiro K, Oishi A, Ooto S, Nakata I, Akagi-Kurashige Y, Miyake M, Yoshimura N, “Macular Choroidal Thickness Measured by Swept Source Optical Coherence Tomography in Eyes with Inferior Posterior Staphyloma.”, *Invest Ophthalmol Vis Sci*. 2012 Oct 25. pii: iovs.12-9952v1.
- “ Ellabban AA, Tsujikawa A, Matsumoto A, Yamashiro K, Oishi A, Ooto S, Nakata I, Akagi-Kurashige Y, Miyake M, Elnahas HS, Radwan TM, Zaky KA, Yoshimura N., “Three-Dimensional Tomographic Features of Dome-Shaped Macula by Swept-Source Optical Coherence Tomography.” *Am J Ophthalmol*. 2012 Nov 3.

- “ Ruiz-Moreno JM, Flores-Moreno I, Lugo F, Ruiz-Medrano J, Montero JA, Akiba M., “Macular choroidal thickness in normal pediatric population measured by Swept-Source Optical Coherence Tomography.”, *Invest Ophthalmol Vis Sci.* 2012 Dec 18.
- “ Ruiz Moreno JM, “Choroidal Imaging With Swept-source Optical Coherence Tomography”, *Retina Today*, November/December 2012
- “ Kyoko Ohno-Matsui, “The Shape of the Sclera Using Swept-source OCT in Eyes With Pathologic Myopia”, *Retinal Physician*, Jan 2013
- “ Dr Florence Coscas, Dr Eric Souied et. al., “Swept Source OCT versus Spectral Domain-EDI-OCT dans les DMLA Revue des outils de suivi”, *Pratiques en Ophtalmologie* • Février 2013 • vol. 7 • numéro 61 [French]
- “ Itakura H, Kishi S, Li D, Akiyama H., “Observation of posterior precortical vitreous pocket using swept-source optical coherence tomography.”, *Invest Ophthalmol Vis Sci.* 2013 May 3;54(5):3102-7
- “ Kaweh Mansouri, Robert N. Weinreb, “ Evaluation of retinal and choroidal thickness by swept source optical coherence tomography: repeatability and assessment of artifacts.”, *AJ Ophthalmol* 2014;157:1022–1032. 2014 by Elsevier
- “ Naoko ueda-arakawa, sotaro ooto,nagahisa yoshimura, “Macular Choroidal Thickness and Volume of Eyes With Reticular Pseudodrusen Using Swept-Source Optical Coherence tomography”, *AJ Ophthalmol* 2014;157:994–1004.
- “ Sergio Copete, José M Ruiz-Moreno, “Direct comparison of spectral-domain and swept-source OCT in the measurement of choroidal thickness in normal eyes”, *Br J Ophthalmol* 2014;98:334–338
- “ Kaweh mansouri, robert n. weinreb, “ Improved visualization of deep ocular structures in glaucoma using high penetration optical coherence tomography”, *Expert Rev. Med. Devices* 10(5), 621–628 (2013)
- “ Kyoung Min Lee & Se Joon Woo & Jeong-Min Hwang, “ Evaluation of congenital excavated optic disc anomalies with spectral-domain and swept-source optical coherence tomography”, *Graefes Arch Clin Exp Ophthalmol* DOI 10.1007/s00417-014-2680-9
- “ Munemitsu Yoshikawa, Tadamichi Akagi, Nagahisa Yoshimura, “Alterations in the Neural and Connective Tissue Components of Glaucomatous Cupping After Glaucoma Surgery Using Swept-Source Optical Coherence Tomography”, *Invest Ophthalmol Vis Sci.* 2014;55:477–484.
- “ Alexandre Pedinielli, Eric H. Souied, Giuseppe Querques,” *In Vivo* Visualization of Perforating Vessels and Focal Scleral Ectasia in Pathological Myopia”, *Invest Ophthalmol Vis Sci.* 2013;54:7637–7643.
- “ Toshihiko Nagasawa, Yoshinori Mitamura, Hitoshi Tabuchi, “ Macular Choroidal Thickness and Volume in Healthy Pediatric Individuals Measured by Swept-Source Optical Coherence Tomography”, *Invest Ophthalmol Vis Sci.* 2013;54:7068–7074
- “ Kohei Takayama, Nagahisa Yoshimura, “ Three-Dimensional Imaging of Lamina Cribrosa Defects in Glaucoma Using Swept-Source Optical Coherence Tomography”, *Invest Ophthalmol Vis Sci.* 2013;54:4798–4807
- “ Yukiko Matsuo, Taiji Sakamoto, “Comparisons of Choroidal Thickness of Normal Eyes Obtained by Two Different Spectral-Domain OCT Instruments and One Swept-Source OCT Instrument”, *Invest Ophthalmol Vis Sci.* 2013;54:7630–7636.
- “ Jorge Ruiz-Medrano, Jos’e M. Ruiz-Moreno, “Macular Choroidal Thickness Profile in a Healthy Population Measured by Swept-Source Optical Coherence Tomography”, *Invest Ophthalmol Vis Sci.* 2014;55:3532–3542
- “ “Danjie Li, Shoji Kishi, Hirotaka Itakura, ““Posterior Precortical Vitreous Pockets and Connecting Channels in Children on Swept-Source Optical Coherence Tomography””, *Invest Ophthalmol Vis Sci.* 2014;55:412–2416.”
- “ Hirotaka Itakura, Shoji Kishi, “Vitreous Changes in High Myopia Observed by Swept-Source Optical Coherence Tomography”, *Invest Ophthalmol Vis Sci.*2014;55:1447–1452.
- “ Kaweh Mansouri, Felipe A. Medeiros, Robert N. Weinreb, “Assessment of Choroidal Thickness and Volume during the Water Drinking Test by Swept-Source Optical Coherence Tomography”, *Ophthalmology* 2013;120:2508-2516
- “ Karen B. Schaal, Claudine E. Pang, Michael Engelbert, “The Premacular Bursa’s Shape Revealed In Vivo by Swept-Source Optical Coherence Tomography”, *Ophthalmology* 2014;121:1020-1028
- “ LS Lim, G Cheung and SY Lee, “Comparison of spectral domain and swept-source optical coherence tomography in pathological myopia”, *Eye* (2014) 28, 488–491
- “ Anna L. Silverman, Andrew J. Tatham, Felipe A. Medeiros, Robert N. Weinreb, “Assessment of Optic Nerve Head Drusen Using Enhanced Depth Imaging and Swept Source Optical Coherence Tomography”, *Silverman et al: J Neuro-Ophthalmol* 2014; 34: 198-205

- “ Kaori Sayanagi, Fumi Gomi, Masahiro Akiba, Miki Sawa, Chikako Hara, Kohji Nishida, “ En-face high-penetration optical coherence tomography imaging in polypoidal choroidal vasculopathy”, doi: 10.1136/bjophthalmol-2013-304658
- “ Zofia Michalewska, Janusz Michalewska, Zofia Nawrocka, Karolina Dulciewska-Cichecka, Jerzy Nawrocki, “Suprachoroidal Layer and suprachoroidal space delineating the outer margin of the choroid in swept-source optical coherence tomography”, *Retina* 0;1-6,2014
- “ Zofia Michalewska, Janusz Michalewska, Ron A. Adelman, EWA Zawlsak, Jerzy Nawrocki, “Choroidal Thickness Measured with Swept Source Optical Coherence Tomography before and after vitrectomy with interal limiting membrane peeling for idiopathic epiretinal membranes”, *Retina* 0;1-5,2014
- “ Janusz Michalewska, Zofia Michalewska, Zofia Nawrocka, Maciej Bednarski, Jerzy Nawrocki, “Correlation of choroidal thickness and volume measurements with axial length and age using swept-source optical coherence tomography and optical low-coherence reflectometry”, *BioMed research international* volume 2014, article ID 639160
- “ Tanawade RG, Muqit MM Mcleod D, Stanga PE. “Swept-source optical coherence tomography imaging in conservative and surgical management of premacular haaemorrhages showing inflammatory response”, *Clin Experiment Ophthalmol.* 2015 Jan;43(1):77-79
- “ Muqit MM, Stanga PE. “Swept-source optical coherence tomography imaging of the cortical vitreous and the vitreoretinal interface in proliferative diabetic retinopathy: assessment of vitreoschisis, neovascularisation and the internal limiting membrane” *Br J Ophthalmol.* 2014 Jul;98(7):994-7
- “ Stanga PE, Sala-Puigdollers A, Caputo S, Jaberansari H, Cien M, Gray J, D’Souza Y, Charles SJ, Biswas S, Henson DB, McLeod D, “In vivo imaging of cortical vitreous using 1050nm swept source deep range imaging optical coherence tomography”, *Am J Ophthalmol.* 2014 Feb,157(2):397-404e2.
- “ Sam Razavi, Eric H. Souied, Edoardo Cavallero, Michel Weber, Giuseppe Querques, “Assessment of Choroidal Topographic Changes by Swept Source Optocal Coherence Tomography After Photodynamic Therapy for Central Serous Chorioretinopathy”, *Am J Ophthalmol* 2014;157:852-860

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Topcon Espanha S.A.