

Specifications

Refractive Power Measurement	Spherical refractive power: -25D to +22D (0.12D/0.25D steps) Cylindrical refractive power: OD to ±10D (0.12D/0.25D steps) Astrigmatism/axial angle: 1° to 180° (in 1° or 5° steps) Minimal measurable pupil diameter: <2 mm
Corneal Curvature Measurement	Corneal curvature radius: 5.00 to 10.00mm (0.1mm steps) Corneal refractive power: 67.50D to 33.75D (0.12D/0.25D steps) (where corneal refractive power = 1.3375) Corneal astigmatic refractive power: OD to ±10D (0.12D/0.25D steps) Corneal astigmatism/axial angle: 1° to 180° (1°/5° steps)
PD Measurement Range	20mm to 85mm (0.5mm steps)
Data Transport Terminal	USB (Import) / RS-232C (Import/Export) / LAN (Export)
Dimensions	321x340mm (W) x 523x59mm (D) x 490-520 mm (H)
Weight	18 kg
Power Supply	100-240V AC, 50-60Hz, 55-80VA

1. Spherical/Cylindrical refractive power = (spherical refractive power + cylindrical refractive power + astigmatism refractive power) × 1.3375



2. Please ask your Dealer or distributor for details and availability of other systems and software.

3. Connection with other software and other systems is not certified by TÜV.



*Subject to change. Refer to the latest specification sheet for detailed information.

**Please refer to the User Instruction Manual for more information on how to obtain the best results with this instrument.

***Please refer to the User Instruction Manual for more information prior to operation.



TOPCON HARMONIC

JACOB

TESTING INSTITUTE

TOPCON TOPCON CORPORATION

TOPCON MEDICAL SYSTEMS, INC.

TOPCON HEALTHCARE SOLUTIONS, INC.

TOPCON CANADA, INC.

TOPCON MEDICAL SYSTEMS, LTD.

TOPCON EUROPE MEDICAL BV.

TOPCON MEDICAL SYSTEMS, LTD.

TOPCON MEDICAL SYSTEMS, LTD.